

PromoChrom Technologies

For automation of sample preparation



Approaches for Automation of sample preparation


Direct automation of manual procedures: Offline SPE




Integration with instrumental analysis: Online SPE



Core technologies: design and manufacture of special valves




A multi functional valve for
multiple channel fluid handling
US patent number 8813785 B2




Status 1. Pumps are connected to
samples.

Valve for multi channel fluid handling




Status 2. Pumps are connected to solvents.




Status 3. Pumps are connected to columns.

An application example of the multifunctional valve



Top: Module for one channel in a typical multi channel SPE and components in a 4 channel SPE

Right: Structure of a 6-channel SPE from PromoChrom. Only two valves are used regardless number of channels




Product 1: SPE-03 8-channel SPE system

1. Process 8 samples in parallel mode
2. Suitable for both large volume water samples and small volume food sample extracts
3. Much smaller size and weight than other multi channel SPE (12 KG)
4. Extra functions: online evaporation after cleanup
5. Robust and easy operation thanks to highly integrated valves



Product 1: SPE-03 8-channel SPE system




The whole system has no XYZ motion components

Product 2: SPE-04 online/offline SPE system

1. Perform both offline SPE and online SPE
2. Use SPE columns available on the market
3. Online derivatization using two reagents and heating after SPE cleanup
4. Easy integration with LC and LC-MS
5. Easy method development




Product 2: SPE-04 online/offline SPE system



Product 3: SPE-04+ 4 channel SPE system

1. Process 40 samples at a speed 4 times of a conventional SPE system
2. Use 16X100 tube for collection
3. Online mix of two solvents for gradient elution.
4. Needle will follow the level of liquid while taking sample to avoid contamination of plunger
5. Simple structure based on the multifunctional valve
6. Small size with built in touch screen computer



Product 3: SPE-04+ 4 channel SPE system

User interface
of the
software

The screenshot shows the main software window titled "SPE-04 plus VER2015.01" by "PromoChrom Technologies Ltd.". It features a menu bar with "Save", "Save as", "Method", "Start", "Pause", and "Setting". Below the menu is a list of methods (01-10) with checkboxes and file names. To the right is a table with columns: Inlet 1, Inlet 2, In2 %, Outlet, Flow, and Volume. At the bottom of the table are buttons for "Ready", "Append", "Insert", "Edit", and "Delete".

| Method | File Name | Inlet 1 | Inlet 2 | In2 % | Outlet | Flow | Volume |
|--|--------------|---------|---------|-------|--------|------|--------|
| <input checked="" type="checkbox"/> 01 | caffeine.txt | | | | | | |
| <input type="checkbox"/> 02 | caffeine.txt | sol 1 | sol 2 | 10 | column | 5 | 5.0 |
| <input checked="" type="checkbox"/> 03 | sugar.txt | sol 1 | sol 2 | 70 | column | 5 | 4.0 |
| <input type="checkbox"/> 04 | sugar.txt | samp | sol 2 | 0 | column | 3 | 6.0 |
| <input type="checkbox"/> 05 | clean.txt | sol 2 | sol 2 | 0 | column | 5 | 5.0 |
| <input type="checkbox"/> 06 | method1.txt | N2 | sol 2 | 0 | column | 5 | 10.0 |
| <input type="checkbox"/> 07 | method1.txt | sol 2 | sol 1 | 40 | frac 1 | 5 | 5.0 |
| <input type="checkbox"/> 08 | method2.txt | sol 2 | sol 1 | 60 | frac 2 | 5 | 5.0 |
| <input type="checkbox"/> 09 | method3.txt | sol 2 | sol 1 | 90 | frac 3 | 5 | 5.0 |
| <input type="checkbox"/> 10 | method3.txt | | | | | | |

Panel for method
parameter entry

The dialog box is titled "Panel for method parameter entry" and contains the following fields and controls:

- Three dropdown menus for "Inlet 1", "Inlet 2", and "Outlet".
- "In2 ratio" field with a value of 1 and up/down arrows.
- "Volume" field with a value of 0 and up/down arrows.
- "Flow rate" field with a value of 0 and up/down arrows.
- "In 0.1 mL" field with a value of 1 and up/down arrows.
- "OK" and "Cancel" buttons at the bottom.

Product 3: SPE-04+ 4 channel SPE system

Table 3. SPE Method (“PAM8 L5.spe”)

| # | Step | Source | Output | Vol. (ml) | Speed (ml/min) | Liquid sense |
|---|---------------|------------------|---------|-----------|----------------|--------------|
| 1 | Condition | 0.5% TEA in MeCN | Waste2 | 5 | 30 | No |
| 2 | Condition | 0.5% TEA in MeCN | Waste2 | 5 | 30 | No |
| 3 | Load | Sample | Fract1 | 5.2 | 9 | No |
| 4 | Collect | MeCN | Fract1 | 5 | 9 | No |
| 5 | Purge-cannula | MeCN | Cannula | 2 | 30 | No |
| 6 | Purge-cannula | MeCN | Cannula | 2 | 30 | No |

Transfer from a Rapidtrace method to a SPE-04+ method.

| Inlet 1 | Inlet 2 | Ratio | Outlet | Flow rat | Volume |
|---------|---------|-------|--------|----------|--------|
| Solv 2 | Solv 2 | 0 | column | 8 | 10 |
| sample | Solv 2 | 0 | Frac 1 | 8 | 5 |
| Solv 1 | Solv 2 | 0 | Frac 1 | 8 | 5 |
| Solv 1 | Solv 2 | 0 | wash | 10 | 10 |



Product 4: LC-03 online SPE system

1. Dedicated to water analysis
2. Achieve better sensitivity and much faster speed than offline SPE for large volume water extraction
3. Sample extraction and LC analysis are performed parallel (no waste time)
4. Computer is not necessary
5. Include SPE column, LC column and methods as a ready solution




Product 4: LC-03 online SPE system

Working principle



Product 4: LC-03 online SPE system

Compare with other solutions on the market



LC-03: The pump can elute the SPE column from two directions. I can easily remove particle without using filter and keep a long lifetime for the SPE column.

The sample loading and elution must be in the same direction. Trapped particles is difficult to remove.

Product 4: LC-03 online SPE system

Example for analysis of benzopyrene in tap water with UV detector

