



UV-Professional

Analysis software

User's Manual

Contents

1. Introduction	1
1.1 Main function.....	1
1.2 Spectra treatment.....	2
2. Installation.....	3
2.1 System configuration	3
2.2 Installation UV Professional	3
2.3 Uninstall UV Professional	4
2.4 Running UV Professional	4
2.5 How to connect the instrument and PC.....	5
2.6 User information	5
2.7 Connect / Release host	5
3. How to use the software	6
3.1 Main interface.....	6
3.2 Menu bar and toolbar	6
Toolbar button list:	6
4. Operation	9
4.1 Photometry	9
4.1.1 Measurement	9
4.1.2 How to test	9
4.2 Quantitative	10
4.2.1 Measurement	10
4.2.2 Coefficient	11
4.2.3 The calibration standard samples	12
4.2.4 Rename ample	13
4.2.5 Save data	13
4.3 Time Scan.....	14
4.3.1 Measurement	14
4.3.2 Change display mode	15
4.3.3 Save data	16
4.4 Wavelength scan.....	17
4.4.1 Measurement	17
4.4.2 Switch display mode.....	19
4.4.3 Display area to enlarge.....	19
4.4.4 Peak search	19
4.4.5 Save data	19
4.4.6 Save data for the picture or text file	20
4.5 Multi-Wavelength	21
4.5.1 Test	21
4.5.2 Rename a sample	22
4.5.3 Save data	23
4.6 DNA / protein analysis	24

4.6.1 Test	24
4.6.2 Renaming sample	25
4.6.3 Save data	25
5. Other features	26
5.1 Set wavelength	26
5.2 Set lamp switch wavelength	26
5.3 Turn on/off tungsten lamp	26
5.4 Turn on/off deuterium lamps.....	26
6 File operations.....	27
6.1 Save the test data.....	27
6.2 Open Test	27
6.3 Print test reports.....	27
Appendix 1	28

1. Introduction

UV Professional has the main function of Photometry, Quantitative analysis, Kinetics, Wavelength Scanning, Multi-wavelength analysis, DNA/Protein analysis. Software attached several additional features, for example: spectrum computing(spectrum addition, spectrum subtraction, spectrum multiplication and spectrum division). spectrum derivative and spectrum compare and so on.

1.1 Main function

Photometry

- Two Test Mode: Absorbancy Data(Abs)/Transition Data(%T) in the custom wavelength.

Quantitative Analysis

- Create standard curve with calibration standard samples ,or with user setup coefficient method.
- Up to 20 standard samples to the standard calibration curve or directly enter the standard curve coefficient.
- 3 methods to standard curve fitting: zero-crossing first-order linear fitting, First-order linear fitting, second-order linear fit).

Time Scan

- Custom scan time interval:0.5 seconds, 1.0 seconds, 2.0 seconds, 5.0 seconds, 10.0 seconds, 30.0 seconds, and 1 Minute.
- 2 display mode: Absorbancy Data (Abs), Transition Data (%T).

Wavelength scanning

- Custom scan wavelength interval:0.1nm, 0.2nm, 0.5nm, 1.0nm, 2.0nm and 5.0nm.
- 3 test mode: Absorbancy Data (Abs), Transition Data (%T) and energy.
- 3 display mode: Absorbancy Data (Abs), Transition Data (%T) and energy.
- system baseline can be stored.

Multi-wavelength analysis

- up to 15 wavelengths

DNA / Protein Analysis

- built-in 2 analytical methods
- Test coefficient can be customized

1.2 Spectra treatment

Spectral data show

- Move the cursor to the spectrum position, the corresponding data will be displayed

Peak valley automatically search

- Automatically search the data peak and valley, and the data will be showed in the peak data list, the corresponding data remark shown in the picture.

Zoom in/out spectrum view

- Change the X Position/Y Position scale data to zoom in or zoom out the spectrum view.

Derivative spectra

- Calculate and display first-order to fourth-order derivative spectrum view. In the absorbance mode, the derivative spectrum view is an extremely effective method.

Spectrum computing

- Spectrum addition, spectrum subtraction, spectrum multiplication, spectrum division and spectrum derivative.

Dark Current Detect

- Retest the dark current data of instrument.

2. Installation


This chapter will show you how install UV Professional to your personal computer.

2.1 System configuration


- Pentium or faster processor
- CD-ROM drive
- 2 USB interface
- 32MB RAM (recommend 256MB or more)
- 50MB or more hard disk space
- Microsoft Windows 2000 or Windows XP operating system

2.2 Installation UV Professional


1. Put UV Professional CD into your CD-ROM drive.
2. In the CD-ROM root directory, double-click Setup.exe to start the installation progress, click <Next> in the setup dialog, then select <accept this agreement...>, continue installation progress.



3. Clicking <change...> button custom the installation directory, click <Setup>, then to start to copy files to your computer.



4. Click <Finished> to complete UV Professional installation process, setup progress automatically install UV key driver and communications USB port driver.




5. The installation process finished.

2.3 Uninstall UV Professional

There are 3 ways to safely uninstall UV Professional

1. In <Control Panel> and <Add or Remove Programs> dialog, , select <UV Professional>, click the <uninstall> button.
2. In the start menu, <program> select <UV Professional> and <uninstall the UV professional>,the uninstall progress will automatically uninstall the UV professional software.



2.4 Running UV Professional

After UV Professional progress is finished, insert UV Key to your PC, connect the instrument and PC with link cable..



There are two ways to run the UV Professional software

1. Double-click on the desktop  [UV Professional] icon.
2. In the [Start] menu-> [All Programs] -> [UV Professional] -> [UV Professional], click [UV Professional] to start run UV Professional.

2.5 How to connect the instrument and PC.

1. Connect the instrument and PC with link cable.
2. click connect button, software will automatically search the communication port to the instrument..

2.6 User information

Select the main menu <View>-><Option>-><Information>, enter the operator and the department information, click OK to save the change, this information will used on print test reports.

2.7 Connect / release host

1. Turn on the instrument power, the instrument will run diagnostic program and warm up procedure. After this, the instrument will be in main menu.


2. Click the shortcut on the toolbar button in the [ Connect/Release] or select the main menu <operation>->[Connect/Release], software will automatically connect with the instrument, the toolbar button display the icon [ Connect/Release] , click this button to release the instrument.

3. How to use the software

This chapter shows you how to use UV Professional

3.1 Main interface

The main interface in the software startup:



3.2 Menu bar and toolbar

Software has the Menu bar and tools bar to provide user the easy ways use this software. At the same time, pop-up menu of right button include most commonly used functions to speed up user operation.

Toolbar button list:

Main	list			introduction
------	------	--	--	--------------

File(F)	New file		[New file]	selection the new file type list, user can choose to photometry, quantitative, time scan, wavelength scan, multi-wavelength, DNA / protein.
			[Photometry]	Create new photometry file.
			[Quantitative]	Create new quantitative file.
			[Time Scan]	Create new time scan file.
			[Wavelength scan]	Create new wavelength scan file.
			[Multi-wavelength]	Create new multi-wavelength file.
			[DNA/protein]	Create New DNA/Protein file.
	Open		[Open...]	Open the test data files.
	Save		[Save...]	Save the test data files.
	Print		[Print...]	Print test report.
Operate	Connect/release	 	[Connect /Release]	Connect or release the instrument.
	Test		[Start test]	start test.
	Stop		[Stop testing]	stop test.
	Tungsten lamp	 	[Turn on/off tungsten lamp]	Turn on/off tungsten lamp.
	Deuterium lamps	 	[Turn on/off deuterium lamp]	Turn on/off deuterium lamps.
	Lamp switch wavelength		[Lamp switching wavelength]	Set lamp switch wavelength.
	Modify		[Modified]	Modify test record.

	Delete		[Remove]	Remove test record
	System baseline		[Create system baseline]	Create the system baseline
	Set blank		[Set blank]	Set blank(0.000Abs/100.0% T)
	Set wavelength		[Set wavelength]	Set wavelength
	Option		[Option...]	setup the option
Spectrum Operation	peak search		[Peak search]	Search the data peak and valley.
	Set Peak value		[Set Peak value]	set peaks and valley value
	Computing spectrum		[Spectrum smoothing]	spectrum smoothing
			[Spectrum add]	the two spectra of addition operations
			[Spectrum subtraction]	two spectral subtraction operator
			[Spectrum multiplied]	two spectral multiplication
			[Spectral division]	the two spectra computing division
			[Spectrum Derivative]	spectra derivative operator
	cascading window		[Cascade]	Cascade window
Windows			[Tile Horizontally]	Tile horizontally window
			[Tile Vertically]	Tile vertically window
Help	help		[Help...]	Help books.
	About		[About...]	Software version information and product ID information.

4. Operation


This chapter show user the software main function: photometry, quantitative, time scan, the wavelength scan, multi-wavelength, DNA/protein.

4.1 Photometry

Photometry have two Test Mode: Absorbancy Data(Abs)/Transition Data(%T) in the custom wavelength.

4.1.1 Measurement

1. Click the shortcut toolbar [Photometry], or <New files> -><Photometry > to create a new photometry file.



4.1.2 How to test

1. Put the sample into the light path, select the menu <operation> <Set Blank> or click the shortcut toolbar buttons [Set Blank],

instrument will be set to the operating wavelength, ans set blank 0.000 Abs/100.0 %T.

2. Put the samples into light path, select the main menu <operation>


<start test> or click the shortcut toolbar buttons [ start test] to measurement the current Absorbancy Data(Abs) and Transition Data(%T), then display the data in the data list.

4.2 Quantitative


Quantitative have two methods: standard curve method and coefficient method to measure the concentration.

4.2.1 Measurement

1. Click the shortcut toolbar [ Quantitative], or <New files>-><Quantitative > to create a new quantitative file.




2. Shortcut toolbar [ Option] to setup quantitative methods, concentration units, the number of wavelengths, coefficient. Click <OK> to save configuration.




4.2.2 Coefficient

1. Select Coefficients methods, third-order fitting, single wavelength 500nm. The following diagram is configuration:



2. Put the sample into the light path, select the menu <operation> <Set Blank> or click the shortcut toolbar buttons [Set Blank], instrument will be set to the operating wavelength, and set blank 0.000 Abs/100.0 %T. Coefficient formula is shown as follows.



3. Put the samples into light path, select the main menu <operation> <start test> or click the shortcut toolbar buttons [start test] to measurement the current absorbance data, calculate the concentration within the above formula, then display the data in the data list.

4.2.3 The calibration standard samples

1. Select standard samples calibration option, change the sample number(maximum number is 20), click <OK> button to complete the setup.
2. In the standard sample data list table, enter the concentration of each sample.
3. Put blank sample into light path, select the menu <operation>-> <Set Blank> or click the shortcut toolbar buttons [Set Blank(0.000 Abs/100.0%T)], Instrument will be set to the wavelength, the data will be set to 0.000 Abs/100.0%T.
4. Put the standard sample into light path, double-click the corresponding sample the absorbency data table in the standard sample data list. Software will read the the absorbance values of standard sample and display in the list table.
5. Just as the above way, measure all standard samples absorbency.
6. After all standard sample's absorbency is measured, the software

will calculate the curve parameter and display the curve line and formula on the curve picture.

7. Put blank sample into light path, select the menu <operation>-> <Set Blank> or click the shortcut toolbar buttons [ Set Blank(0.000 Abs/100.0%T)], Instrument will be set to the wavelength, the data will be set to 0.000 Abs/100.0%T.

8. Put the samples into light path, select the menu <operation> <start testing> or click the shortcut toolbar buttons [ start test] to measure the sample's absorbency data and calculate the sample concentration, then display the data in the data list.

4.2.4 Rename ample

In the data list to select a sample named Double-click the name to enter the edit state of the sample, enter the name after the carriage return, and then click Save, the sample name on the sample data stored in the same data in a.

4.2.5 Save data


1. Select the menu <File> <Save...> or click shortcut toolbar [ Save...].
2. In the file save dialog, change the file directory and file name, click <OK> button.
3. The quantitative analysis data file is saved with *.qua file extension.

4.3 Time Scan


This chapter show user how to do time scan measurement with the absorbance or transmittance mode.

4.3.1 Measurement

1. Click the shortcut toolbar [ scan time], or <File>-> <time scan>
Create a new time scan file.




2. Click shortcut toolbar [], setup time scan parameters, select display mode, wavelength, start time, end time, time interval.



3. Click <OK> button to save configuration.
4. Put blank sample into light path, select the menu <operation>-> <Set Blank> or click the shortcut toolbar buttons [Set Blank(0.000 Abs/100.0%T)], Instrument will be set to the wavelength, the data will be set to 0.000 Abs/100.0%T.
5. put the sample into light path, select the menu <operation> <start test> or click the shortcut toolbar buttons [start test] to start the time scan measurement.
6. All data will be drawn in the picture and display in the data list.

4.3.2 Change display mode

Shortcut Toolbar [Option], in the option window, select the absorbance mode or transmittance mode.



4.3.3 Save data


1. Select the menu <File> <Save...> or click shortcut toolbar [Save...].
2. In the save file dialog, change the file directory and file name, click <OK> button.
3. Time scan data file Data is saved with an *.tis file extension.

4.4 wavelength scan


This chapter shows user how to do wavelength scan with absorbency, transmittance or energy mode.

4.4.1 Measurement


1. Click on the shortcut toolbar [ wavelength scan] or [file] [wavelength scan] to create new wavelength file.



2. Select the main menu <operation> <set> or click the toolbar shortcut [ Settings], open the Preferences window wavelength scanning.



3. Select display mode (elected to the energy scan mode, the proposal is set to a fixed gain), scanning the wavelength range, coordinates the upper limit, lower limit coordinates and scanning interval.
4. Click OK to complete and exit setup.
5. Will be placed in the reference light path.
6. Click on the shortcut toolbar [ system set up baseline] or <operator> "Create system baseline> Open dialog system baseline to begin scanning system baseline or to open the system had previously stored baseline for testing. Click [Start] button to start the baseline correction system. Baseline correction system will take a few minutes, it is proposed that a longer time interval after baseline correction for a system to ensure testing accuracy. If you want to save the current system baseline, click [Save] button to save the current system baseline to a file. Click [Open] button to open the previously stored system baseline file.



7. put reference sample into light path, select the main menu <operation>-><zero / full-scale> or click the shortcut toolbar

buttons [zero / full-scale], software will set blank(0.000Abs/100.0%T) within measurement wavelength range.

8. put test sample into light path, select the main menu <operation>-><start test> or click the shortcut toolbar buttons

[test] to start the test.

4.4.2 Switch display mode

Click shortcut toolbar click [Settings], open the option dialog to choose display mode (transition or absorbance).

4.4.3 Display area to enlarge

In the wavelength scan option dialog, change the top limit and bottom limit, zoom in/out the display spectrum.

4.4.4 Peak search

Click the toolbar shortcut [search peak], after the search finished, the result will be in the list, peak value can be changed.

4.4.5 Save data

1. Select the main menu <File>-><Save> or click shortcut toolbar [Save].

2. In the file save window, select file directory and enter file name, click <OK> button.
3. Data will be saved in the wavelength scan file with file extension of wls wavelength scan data files.

4.4.6 Save data for the picture or text file


Select <File>-><output> to save data for BMP bitmap file or text file txt.

4.5 Multi-Wavelength

This chapter show user how to do multi-wavelengths (up to 15) measurement and to setup the multi-wavelength configuration.

4.5.1 Test

1. Select the main menu <File> <multi-wavelength> or click the toolbar shortcut [ multi-wavelength], create a multi-wavelength measurement.



2. Select the main menu <operation> <option> or click the toolbar shortcut [ Settings], open multi-wavelength analysis of parameter settings window.



3. enter wavelengths number, and enter measurement wavelengths and coefficient:

Results = Abs(wavelength 1)×factor 1 + Abs(wavelength 2)×factor 2 + Abs(wavelength 3)×factor 3 +.....

4. Click <OK> button to complete configuration.

5. Put reference sample into light path, select the main menu <operation> <set zero / full-scale> or click the shortcut toolbar

buttons [zero / full-scale], software will set 0.000Abs/100.0%T for each wavelength. This function will be need several seconds or a few minutes.

6. Put sample into light path, select the main menu <operation> <start test> or click the shortcut toolbar buttons [test] to start test.

4.5.2 Rename a sample

1. select a named text box in the data list table.

2. Double-click, enter the sample name.

4.5.3 Save data


1. Select the main menu <File>-><Save> shortcut toolbar or click  [Save...].
2. In the file save dialog, select file directory and enter file name, click <OK> button.
3. Test Data will be saved as the data file with mul extension.

4.6 DNA / protein analysis


This chapter show user how to do DNA / protein measurement.

4.6.1 Test

1. Select the main menu <File>-><DNA/Protein> or click the toolbar shortcut [ DNA / protein], create a DNA / protein measurement file.



2. Select the main menu <operation>-><Option> or click the toolbar shortcut [ Settings], open the DNA / protein configuration dialog.
3. Select measurement method, the software have two kinds of built-in methods, you can also choose Custom method to enter the test wavelength and coefficient.



4. Click <OK> button to complete configuration.
6. Put reference sample into light path, select the main menu <operation>-><zero / full-scale> or click the shortcut toolbar button [国旗 zero / full-scale], software will set 0.000Abs/100.0%T for each wavelength.
8. put sample into light path, select the main menu <operation> <start test> or click the shortcut toolbar button [播放 test] to start test.

4.6.2 Renaming sample

1. select a named text box in the data list table.
2. Double-click, enter the sample name.

4.6.3 Save data

1. Select the main menu <File>-><Save> shortcut toolbar or click [文件 Save...].
2. In the file save dialog, select file directory and enter file name, click <OK> button.
3. Test Data will be saved as the data file with <.dna> extension.

5. Other features

This chapter shows the other features of UV Professional: set wavelength, set light switch wavelength, turn on/off tungsten lamp, turn on/off deuterium lamp, reset dark current.

5.1 set wavelength

Select main menu <operation>-><set wavelength> or click the toolbar shortcut [ set wavelength], to set work wavelength, click <zero> button to set 0.000Abs/100.0%T for current wavelength.

5.2 Set lamp switch wavelength

Select main menu <operation>-><set lamp switch wavelength> or click the toolbar shortcut [ set lamp switch wavelength] to set set lamp switch wavelength, the default value is between 300nm-400nm.

5.3 Turn on/off tungsten lamp

Select main menu <operation> <Turn on /off tungsten lamp> or click the toolbar shortcut / [ /  Turn on/off tungsten lamp] to turn on or turn off tungsten lamp

5.4 Turn on/off deuterium lamps

Select main menu <operation> <Turn on/off deuterium lamp> or click the toolbar shortcut / [ /  Turn on/off deuterium lamp] to turn on or turn off tungsten lamp

6 File operations

Software file format description

Photometry: *.bas

Quantitative: *. Qua

Time Scan: *. Kin

Wavelength scanning: *. Wls

Multi-wavelength: *. Mul

DNA / protein: *. Dna

System baseline: *. Sbl

6.1 Save the test data

Select main menu <File>-><Save...> or click shortcut toolbar button [Save ...], to display open file save dialog, enter the file name, click <OK> button to save the test data file.

6.2 Open Test

Select main menu <File>-><open> or click the toolbar shortcut button [open ...], to display open file dialog, select the file name, click <OK> button to open the file.

6.3 Print test reports

Select main menu <File> <print> or click the toolbar shortcut [Print ...] to print the test data and spectrum to current printer.

Appendix 1

DNA / Protein measurement methods:

$$C_{DNA} = (Abs_{WL.1} - Abs_{Ref. WL.}) \times F1 - (Abs_{WL.2} - Abs_{Ref. WL.}) \times F2$$
$$C_{Protein} = (Abs_{WL.2} - Abs_{Ref. WL.}) \times F3 - (Abs_{WL.1} - Abs_{Ref. WL.}) \times F4$$

$$\text{Ratio} = (Abs_{[WL.1]} - Abs_{[ref. WL.]}) / (Abs_{[WL.2]} - Abs_{[ref. WL.]})$$

Method 1:

WL. 1 = 260nm, WL. 2 = 280nm, ref. WL. = 320nm

F1=62.9 F2=36.0 F3=1552 F4=757.3

$$C_{DNA} = (Abs_{[260nm]} - Abs_{[320nm]}) \times 62.9 - (Abs_{[280nm]} - Abs_{[320nm]}) \times 36.0$$

$$C_{Protein} = (Abs_{[280nm]} - Abs_{[320nm]}) \times 1552 - (Abs_{[260nm]} - Abs_{[320nm]}) \times 757.3$$

$$\text{Ratio} = (Abs_{[260nm]} - Abs_{[320nm]}) / (Abs_{[280nm]} - Abs_{[320nm]})$$

Method Two:

WL. 1 = 260nm, WL. 2 = 230nm, ref. WL. = 320nm

F1=49.1 F2=3.48 F3=183 F4=75.8

$$C_{DNA} = (Abs_{[260nm]} - Abs_{[320nm]}) \times 49.1 - (Abs_{[230nm]} - Abs_{[320nm]}) \times 3.48$$

$$C_{Protein} = (Abs_{[230nm]} - Abs_{[320nm]}) \times 183 - (Abs_{[260nm]} - Abs_{[320nm]}) \times 75.8$$

$$\text{Ratio} = (Abs_{[260nm]} - Abs_{[320nm]}) / (Abs_{[230nm]} - Abs_{[320nm]})$$

AOE INSTRUMENTS (SHANGHAI) CO., LTD All Rights Reserved