

# UV-Win GLP

## SOFTWARE

**UV-Win is a powerful, intuitive Software product used for connectivity to the PG Instruments range of bench top UV-Vis Spectrophotometers.**



**UV-Win GLP offers all of the features and functionality of UV-Win whilst also offering extensive administrative capabilities along with a detailed audit trail.**

The UV-Win software offers complete instrument control along with data acquisition and a whole host of mathematical tools for interpretation of measurement results. The UV-Win software is separated into four key workspaces:

- Spectral Analysis.
- Quantitative Analysis.
- Kinetic Analysis.
- Photometric Analysis.


### SPECTRUM WORKSPACE

- Use the spectrum workspace to scan across a userdefined spectral range measuring in either absorbance or transmission.
- Use the "Peak Pick" tool to determine the wavelength at which peaks and valleys have occurred whilst also being able to determine their amplitude.
- View spectral overlay in the 3D display mode.
- Perform 1st, 2nd, 3rd and 4th order differentiation on sample scans for Derivative Spectroscopy.
- Export measurement data into Word, Excel, CSV and ASCII formats.
- Create method files for routine analysis whilst also being able to save measurement data.



### QUANTITATIVE WORKSPACE

- Use the Quantitative workspace to determine the concentration of unknown samples.
- Create a calibration curve using a series of standard solution or by directly entering the coefficients for the calibration curve equation.
- Specify 1st, 2nd, 3rd and 4th order correlation for the best calibration curve fit.
- Set Quality Control monitors to take user specified action in the event of measurement results falling outside user defined limits.
- Export measurement data into Word, Excel, CSV and ASCII formats.
- Create method files for routine analysis whilst also being able to save measurement data.



Set quantitation parameters

Accessories:  Quality control  Simple Calculation

Measurement:  Calibration curve  Instrument

Curve equation:  $Abs = K1 * Conc + K0$  Equation Order: 1st

Input Type:  Conc  Values

Conc. unit: mg/l  Blank Extension CheckBox  Zero Interceptor

Curve evaluation: None

Calibration method:  Concentration method  Coefficient method


K1: [ ] K0: [ ] K2: [ ]

K3: [ ] K4: [ ]

Buttons: OK, Cancel, Help


## KINETIC WORKSPACE

- Monitor the change of Absorbance or Transmission as a function of time for Enzyme type reactions.
- Use in conjunction with a Flowcell for sample introduction or Peltier water circulator for temperature control.
- Specify data intervals and acquisition time for up to 24 hour reactions.
- Export measurement data into Word, Excel, CSV and ASCII formats.
- Create method files for routine analysis whilst also being able to save measurement data.



## PHOTOMETRIC WORKSPACE

- Perform a series of sequential fixed wavelength measurements in either Absorbance or Transmission.
- Automate sample measurements by configuring the instrument cell changer.
- Calculate concentration of unknown samples quickly using the "Simple Calculation" option where complete calibration is not required.
- Automatically calculate statistics like standard deviation, relative standard deviation, and averages.
- Export measurement data into Word, Excel, CSV and ASCII formats.
- Create method files for routine analysis whilst also being able to save measurement data.



No.	Wavelength	Abs	Trans	Conc	Std Dev	Rel Std Dev
1	1.0	0.000	0.999	0.000	0.000	0.000
2	1.1	0.000	0.999	0.000	0.000	0.000
3	1.2	0.000	0.999	0.000	0.000	0.000
4	1.3	0.000	0.999	0.000	0.000	0.000
5	1.4	0.000	0.999	0.000	0.000	0.000
6	1.5	0.000	0.999	0.000	0.000	0.000
7	1.6	0.000	0.999	0.000	0.000	0.000
8	1.7	0.000	0.999	0.000	0.000	0.000
9	1.8	0.000	0.999	0.000	0.000	0.000
10	1.9	0.000	0.999	0.000	0.000	0.000



## REPORTING

- Produce reports for photometric, spectrum, kinetic and quantitative measurement data.
- Include or remove spectra, calibration curves along with samples measurement tables.

