

Multi-Tube Vortex Mixer

HD2500



Features

HD2500 can process 50 samples at once. Several accessories can be selected to meet the mixing needs ➤

LCD display, micro-processor control

- Equipped with 12mm a foam tube rack and tray mats
- The panel is simple and the microprocessor precisely controls the time and speed
- The maximum speed can reach 2500 rpm, well mixing effect

Specifications

Speed	500 ~ 2500rpm
Speed Accuracy	±1rpm
Amplitude	3.6mm
Timer	0 ~ 99m59s
Max. Load Capacity	4.5Kg
Platform Dimension (mm)	310(w) x 180(D)
Power	75W
Voltage	AC110~230V; 50~60Hz
Weight	14.7kg
Dimension (mm)	426(W) x 246(D) x 474(H)

Accessories

	HD2500/X1	Attachment-X1	50 x Ø12mm
	HD2500/X2	Attachment-X2	50 x Ø13mm
	HD2500/X3	Attachment-X3	50 x Ø16mm (15ml Tube)
	HD2500/X4	Attachment-X4	15 x Ø25mm
	HD2500/X5	Attachment-X5	15 x Ø29mm (50ml Tube)
	HD2500/X6	Attachment-X6	Mat 305 x 178.5 x 25 (mm)