

# VACUUM DRYING OVEN

## LAS-VD SERIES

LAS-VD series vacuum drying ovens are designed for drying heat-sensitive, easily decomposed and easily oxidized substances through perfect heat conduction. When working, it can maintain a certain degree of vacuum in the working room, and can fill the interior with inert gas, and can also quickly dry items with complex components.



## PRODUCT FEATURES

- PID temperature control technology, accurate and reliable temperature control, LCD screen displays the set temperature and current temperature, clear and intuitive, with timing, over-temperature alarm, automatic error correction and other functions;
- The partition can be loaded and unloaded freely, easy to clean;
- The closed tightness of the box door can be adjusted, and the integrally formed silicone sealing strip ensures a high vacuum in the box, and inert gas can be filled into the working room (the inflation pressure should not exceed 0.1MPa);

Product Model	LAS-VD-20	LAS-VD-50	LAS-VD-90	LAS-VD-210
Image				
Voltage (V)	AC220			
Frequency (Hz)	50			
Power (KW)	0.6	1	1.8	3
Display Mode	LCD			
Vacuum (pa)	133			
Temperature Control Range (°C)	RT+10-280			
Temperature Control Accuracy (°C)	±1			
Timing Range (min)	0-9999			
Working Size (mm)	300*300*275	370*415*345	450*450*450	600*550*650
Effective Volume (L)	22	52	91	214
Standard Number of Shelves	2		3	
Bearing of Shelf (kg)	15			